

Table S1. Comparison of thyroid function test (TFT) results obtained with the Siemens ADVIA XP testing platform before and after polyethylene glycol (PEG) treatment

	FT3 [pmol/L]	FT4 [pmol/L]	TSH [mIU/L]
Original serum	0.49	92.93	> 150
PEG-treated serum	0.47	6.84	> 150
Reference interval	3.5–6.5	11.5–23.5	0.3–5.0

FT3 — free triiodothyronine; FT4 — free thyroxine; TSH — thyroid-stimulating hormone